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
Catalogue No.: V606101-ND V606102-ND V606103-ND	Medical device variants: Infusion set IS-101 Infusion set IS-102 Infusion set IS-103	Catalogue No.: V606103-01ND V606103-03ND V606103-04ND	Medical device variants: Infusion set IS-103 – pack. 4x50 pcs. Infusion set IS-103/K Infusion set IS-103/K – pack. 4x50 pcs
Product type:	medical device		
Design:	sterile, ethylene oxide sterilization		
Target patient:	population of pediatric or adult patients		
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.		
Materials used:	ABS, PE, PVC, PS, PP, PA, RUBBER		
Compatibility:	With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With plastic infusion bags, plastic or glass bottles. With infusion pharmaceuticals without any special requirements for filtration. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz		
Incompatibility:	With infusion pharmaceuticals incompatible with PVC (cytostatics) and with infusion pharmaceuticals sensitive to light.		
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.		
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.		
Warning:	Do not use for blood transfusion. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Infusion bottle of contact with the infusion set until immediately before administration. During the operation, check the patient's condition and functionality sets. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.		
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.		

Intended use

The set is intended for therapeutic application of infusion solutions from bottles or bags using gravitation.

Instruction for use

- Under aseptic conditions, remove the set from the packing.
- Close the hose regulation clamp.
- Take off the cover of the piercing spike and the cover of the cone.
- Inject the piercing spike upright into the plug of the bottle or bag. Hang the bottle or bag onto the stand.
- Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.
- Close the regulation clamp.
- Connect the set on the patient, open the regulation clamp and start the infusion.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Gravity
	The number of drops per milliliter		Filter of liquid with pore size		Medical Device		Date of Manufacture

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

DE
Infusionsset IS-101 / Infusionsset IS-102 / Infusionsset IS-103®
Infusionsset IS-103® - Pack. 4x50 St / Infusionsset IS-103®/K / Infusionsset IS-103®/K - Pack. 4x50 St

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